1 The Fluxgate Magetometer

Chapter 13: Understanding the Fluxgate Magnetometer

Notes

- Teachspin Fluxgate Magnetometer Module:
 - Fluxgate Sensor:
 - Double Solenoid: Two electrically separate wires are wound and stored in a tube for the fluxgate sensor to fit in. This is so we can have up to two independently added magnetic fields.
 - * Solenoid pitch = 2.54 mm, and the expected magnetic field inside the Solenoid is given by

$$B_{\rm ext} = 2\mu_0 ni$$

where i is the current through a solenoid, 2 accounts for the doubled layer, and the turn density is

$$n = 1 \text{ turn}/(0.00254 \text{ m} = 394 \text{ m}^{-1})$$

From this, the magnitude of the external magnetic field given a current i is


$$\frac{B_{\rm ext}}{i} = 2 \cdot (4\pi 10^{-7} \, {\rm Tm/A}) 394 \, {\rm m} = 990 \, \mu {\rm T/A}$$


- Modeling Output: Calibrating the Solenoid
 - * For simple models the 2f-component of the field is linear $A = SB_{\text{ext}}$, but we are actually measuring the magnitude of the spectral component

$$M = S|B_{\text{ext}}|$$

* The geometry of the vector components require us to take a phaser sum which gives us the square magnitude of the measured field

$$M^{s} = S^{2} \left[\left(B_{\text{ext}} + \frac{A_{para}}{S} \right)^{2} + \left(\frac{A_{perp}}{S} \right)^{2} \right] = S^{2} \left[(B_{ext} + a)^{2} + b^{2} \right]$$





- (a) Fluxgate Magnetometer sensor components
- (b) Fluxgate sensor and double solenoid setup

Figure 1.1: (a) Fluxgate Magnetometer sensor components and (b) Fluxgate sensor and double solenoid setup

Experimental Setup

- SR770 Config:
 - FREQ: SPAN 3.125 kHz for observation, 12.2 Hz for measurement; Center Freq. 2 kHZ (second harmonic)
 - MEASure: PSD, Flattop Window
 - Average: 16, ExponentialSOURCE OUT: 1 kHz 1 V

Procedure

- 770: 1 kHz 1 V Sine SOURCE OUT \rightarrow Power Audio Amp. module and adjust gain to 6 V (monitor output with splitter to scope)
- ullet Power Audio Amp. output \to Fluxgate Primary (AC drive) coil
- Secondary (pick-up) $coil \rightarrow SR770$ input
- Place the Fluxgate sensor in the double solenoid pointing straight down
- \bullet Measure second harmonic 2f component of the frequency spectrum
- Add 2.5 V DC current (from 33500 B) to Solenoid A and measure the changes in the 2f component
 - $-\,$ NOTE: Use 50 Ohm Terminator so the measured outut voltage is exactly the displayed voltage due to the internal 50 Ohm impedence of the 33500B
- 770 config change: FREQ SPAN 12.2 Hz
- Solenoid A: Replace DC current with AC ($f_A = 5$ Hz 1 V Sine wave from 33500B) and find the $2f \pm f_A$ components
- Solenoid B: Add a second AC drive ($f_B=2.5~{\rm Hz}~0.5~{\rm V}$ Sine wave from second 33500 B) and measure the four side band frequencies $2f\pm f_A, 2f\pm f_B$

Observations